


Cellometer® Cell Counting Chamber

Cellometer® counting slide consists of two independent enclosed chambers with precisely controlled height. Cell suspension of 15 to 20 μ l is pipette into the chamber through one of the sample introduction ports. Cell counting occurs in the central location of the counting chamber. Precise and complete loading of the counting chamber is not necessary as long as the counting area is covered with cell solution.



CHT4-PD100

Disposable, all plastic cell counting chambers for Cellometer cell counters.
Packed in microscope slide boxes, ready to use.



CHT4-SD100

Disposable, all plastic cell counting chambers for Cellometer cell counters.
Packed with protective film.
Remove protective film before use.

Cat.#	Detailed Description	US Price*	US Price / count
CHT4-PD100-002	Case of 100 slides for 200 counts.	\$300.00	\$1.50
CHT4-PD100-003	Case of 500 slides for 1000 counts.	\$1,350.00	\$1.35
CHT4-PD100-303	Case of 1500 slides for 3000 counts.	\$3,810.00	\$1.27
CHT4-PD100-503	Case of 2500 slides for 5000 counts.	\$5,840.00	\$1.17
CHT4-SD100-004	Case of 750 slides for 1500 counts.	\$1,350.00	\$0.90
CHT4-SD100-304	Case of 2250 slides for 4500 counts.	\$3,810.00	\$0.85
CHT4-SD100-504	Case of 3750 slides for 7500 counts.	\$5,840.00	\$0.78

* Plus shipping and handling.